



NSLIDEG	PNEUMOSLIDE IgG	20 Tests	CE 0318
NSLIDEM	PNEUMOSLIDE IgM	10 Tests	CE 0318
NBACT	PNEUMOBACT	20 Tests	CE 0318

pneumoslides pneumobact

The IFA solution to the serological diagnosis of pneumonia



pneumoslides

Simultaneous detection of the bacteria and viruses most frequently associated with human respiratory infections, using a single slide.

Ref. NSLIDEG

Pneumoslides IgG can be used to carry out a presumptive diagnosis by means of IgG antibody detection. It allows to perform 20 tests.

Ref. NSLIDEM

Pneumoslides IgM allows presumptive diagnosis by measuring IgM responses. It allows to perform 10 tests.

pneumobact

Simultaneous detection of the bacteria most frequently associated with human respiratory infections, using a single slide.

Ref. NBACT

Pneumobact is aimed at the performance of presumptive diagnosis through IgG and IgM antibody detection. It allows to perform 20 test



pneumoslide



VIRCELL
microbiologists

kit contents:

Pneumoslide IgG

20 slides of 10 wells with the following antigens:

- L. pneumophila serogroup 1
- M. pneumoniae
- C. burnetii
- C. pneumoniae
- Adenovirus
- Respiratory syncytial virus
- Influenza A
- Influenza B
- Parainfluenza serotypes 1, 2 and 3
- Cell control

1 vial of PBS.

1 ml of positive control serum.

1 ml of negative control serum.

4 vials of anti-human IgG FITC conjugate.

Mounting medium.

Pneumoslide IgM

10 slides of 10 wells with the following antigens:

- L. pneumophila serogroup 1
- M. pneumoniae
- C. burnetii
- C. pneumoniae
- Adenovirus
- Respiratory syncytial virus
- Influenza A
- Influenza B
- Parainfluenza serotypes 1, 2 and 3
- Cell control

1 vial of PBS.


500 µl of positive control serum.

500 µl of negative control serum.

2 vials of anti-human IgM FITC conjugate.

Mounting medium.

1 vial of anti-human IgG globulin (sorbent).



pneumobact

kit contents:

Pneumobact

10 slides with two rows of 5 wells with the following antigens:

- L. pneumophila serogroup 1
- M. pneumoniae
- C. burnetii
- C. pneumoniae
- C. psittaci

1 vial of PBS.

300 µl of positive control serum for IgG.

300 µl of positive control serum for IgM.


600 µl of negative control serum.

2 vials of anti-human IgG FITC conjugate.

2 vials of anti-human IgM FITC conjugate.

Mounting medium.

1 vial of anti human IgG globulin (sorbent).



Distributor's sticker

NSLIDEG	PNEUMOSLIDE IgG	20 Tests	0318
NSLIDEM	PNEUMOSLIDE IgM	10 Tests	0318
NBACT	PNEUMOBACT	20 Tests	0318