


IDENTIFICATION FROM CULTURE
IN JUST 15 MINUTES!!

VIRAPID LEGIONELLA CULTURE

Immunochromatographic test for the qualitative detection and identification of Legionella genus, *Legionella pneumophila* serogroup 1 and *Legionella pneumophila* serogroups 1-14 in samples from bacterial cultures.

Ref. VR002

25 tests



PRODUCT FEATURES

- Multitest detection and identification of the 3 main groups with clinical interest (LPN1, LPN1-14 and L. genus)
- Rapid detection- Visual reading in just 15 minutes
- User friendly protocol with little manipulation
- All the necessary reagents included in the kit
- Product control line on each device
- Based on specific and highly purified MAb (conjugates and lines)
- 99% sensitivity and 100% specificity average performance*
- Valid for cultures from environmental and clinical samples
- CE marked (IVD)

* Results presented at LEGIONELLA 2009 (7th International Conference. Paris, France)
LPN1 (97% sensitivity, 100% specificity); LPN1-14 (100% sensitivity, 100% specificity);
L. genus (100% sensitivity and 100% specificity). 204 samples tested.